





Sharing raw sequence data to ENA





Data storage & Access



National Food Institute, Technical University of Denmark






Data storage & Access

European Bioinformatics Institute (EBI)








www.ebi.ac.uk/ena/





EMBL-EBI to be HTTPS by default from 1st October

On the 1st October the majority of services hosted on www.ebi.ac.uk will be served over HTTPS by default. Services that are becoming HTTPS by default will automatically redirect users accessing the site on insecure HTTP URLs to secure HTTPS URLs.

Users of EMBL-EBI services may wish to update links, bookmarks or API clients to use the HTTPS URLs.

Search results for **PRJEB22091**

Showing 15 results out of 129 in All results

Filter your results

Source

All results (129)

Nucleotide sequences (129)

Nucleotide sequences (129 results) PRJEB22091 Related data -Danmap 2016 Source: Study ID: PRJEB22091



Navigation	Read Files	Portal	Attributes							
Bulk Download Files 🛕 (Please use Firefox to launch the bulk downloader app.)										
Download: 1 - 127 of 127 results in TEXT										
Select columns										

Showing results 1 - 10 of 127 results

Study accession	Sample accession	Secondary sample accession	Experiment accession	Run accession	Tax ID	Scientific name	Instrument model	Library layout	FASTQ files (FTP)	FASTQ files (Galaxy)	Submitted files (FTP)	Submitted files (Galaxy)	NCBI SRA file (FTP)	NCBI SRA file (Galaxy)	CRAM Index files (FTP)	CRAM Index files (Galaxy)
PRJEB22091	SAMEA104205896	ERS1864914	ERX2148098	ERR2091304	562	Escherichia coli	Illumina MiSeq	PAIRED	File 1	File 1	Fastq file 1 Fastq file 2	Fastq file 1 Fastq file 2	File 1	File 1		
PRJEB22091	SAMEA104205897	ERS1864915	ERX2148099	ERR2091305	562	Escherichia coli	Illumina MiSeq	PAIRED	File 1	File 1 File 2	Fastq file 1 Fastq file 2	Fastq file 1 Fastq file 2	File 1	File 1		
PRJEB22091	SAMEA104205898	ERS1864916	ERX2148100	ERR2091306	562	Escherichia coli	Illumina MiSeq	PAIRED	File 1 File 2	File 1	Fastq file 1 Fastq file 2	Fastq file 1 Fastq file 2	File 1	File 1		
PRJEB22091	SAMEA104205899	ERS1864917	ERX2148101	ERR2091307	562	Escherichia coli	Illumina MiSeq	PAIRED	File 1 File 2	File 1	Fastq file 1 Fastq file 2	Fastq file 1 Fastq file 2	File 1	File 1		
PRJEB22091	SAMEA104205900	ERS1864918	ERX2148102	ERR2091308	562	Escherichia coli	Illumina MiSeq	PAIRED	File 1	File 1	Fastq file 1 Fastq file 2	Fastq file 1 Fastq file 2	File 1	File 1		
PRJEB22091	SAMEA104205901	ERS1864919	ERX2148103	ERR2091309	562	Escherichia coli	Illumina MiSeq	PAIRED	File 1 File 2	File 1		Fastq file 1 Fastq file 2	File 1	File 1		
PRJEB22091	SAMEA104205902	ERS1864920	ERX2148104	ERR2091310	562	Escherichia coli	Illumina MiSeq	PAIRED	File 1	File 1	Fastq file 1 Fastq file 2	Fastq file 1 Fastq file 2	File 1	File 1		
PRJEB22091	SAMEA104205903	ERS1864921	ERX2148105	ERR2091311	562	Escherichia coli	Illumina MiSeq	PAIRED	File 1	File 1	Fastq file 1 Fastq file 2		File 1	File 1		





Submitting data in ENA


- Manual and automated
- Use Webin (interactive web submission system) for new sequencing projects, assembled sequences and annotation
- For other types of data, different channels are available (e.g. datasubs@ebi.ac.uk, FTP, RESTful web-based service)
- www.ebi.ac.uk/ena/about/submit_and_update





Submitting data in ENA


- 1) Install a program to connect to FTP
- https://filezilla-project.org








- 2) Connect to FTP
- Host: webin.ebi.ac.uk
- Username: Your own Webin ID
- Password: Your own password for your ENA account. Then click 'Connect'





• 3) Transfer the genome files to ENA FTP









Submitting data in ENA

- 4) log in ENA with your own webin account
- https://www.ebi.ac.uk/ena/submit/sra/#home













• 5) click 'New Submission' and add 'Center Name'. Then click 'Next'








• 6) Click 'Create a new study'








• 7) Set release date, provide title, and abstract






• 8) 'Select Checklist'








9) Select 'Pathogens Checklists' and then select 'ENA GMI_MDK:1.1'








 10) You have an option to download metadata as excel template or you just click 'Next' to fill out metadata via website








• 11) Add 'title', search 'Escherichia coli' (you have to wait until the system give you organism name to choose), add 'collected_by', add 'collection date'








• 12) Add 'country', 'latitude', 'longtitude', 'environmental_sample', host health state', 'host scientific name', 'host associated' and 'isolate'








13) add number of sample. In this case, it is only 1 sample to add (click '+Add')








 14) Add 'Unique Name' for the submitted sample. In this case, try to add 'Ecoli_001'








• 15) Click type of submitted files. In this case is 'Two Fastq files (Paired)'








16) Add required metadata (*)










• 17) MD5 checksum is also required







• 18) The last step is to click 'Submit'. If something goes wrong, the system will give you error message.

